

PROSAT®

Gowning Room Wipes

Lightweight cellulose/polyester nonwoven wipes

The Gowning Room Wipes are made using a strong cellulose polyester nonwoven fabric that is folded for easy pickup. It is very low in extractable residue. Saturated with a blend of 6% isopropyl alcohol and 94% deionized water these wipes are ideal for glove wipe down during gowning procedures. The wipe is folded for easy pickup.

Packaged in refill packs and a dispensing tub to keep the wipes clean and saturated, the wipes are always conveniently at hand.

These wipes meet the requirements of USP<797> and IEST-CC-RP004.4 for “non-shedding, low-lint, lint-free wipes”.



Features	Benefits
Folded	• For easy pickup
Blended nonwoven fabric	• Very low in extractable residue
Packaged in refill packs and a dispensing tub	• Wipes remain clean and saturated • Wipes are always conveniently at hand

Part No.	Description	Size	Packaging
PSCL0001	PROSAT Gowning Room Wipes, Presaturated with 6% IPA/94% DI Water	5" x 8" (130 x 203 mm)	225/bag; 9 pouches/case

Product Information

Material	54% cellulose/46% polyester
Construction	Hydroentangled staple fibers
Packaging materials	Outer bags (OB1), low density polyethylene (LDPE)  Case (CS), corrugated fiberboard (PAP) 



Technical Data

Attribute (units)	Typical Value	Test Method
Basis weight, nominal; (g/m ²)	53	Contec Method
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2
In deionized water; (g/m ²)	0.04	
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	11.9	
Chloride; (ppm)	18.2	
Particles, readily releasable		IEST-RP-CC004.2, Sec. 5.1
Particles $\geq 0.5\mu\text{m}$; ($\times 10^6/\text{m}^2$)	18.6	
Fibers $\geq 100\mu\text{m}$; ($\times 10^3/\text{m}^2$)	21.7	

Recycle Symbols

PET 

HDPE 

LDPE 

PP 

PAP 

Packaging	EA/PCH	PCH/OB1	EA/CS
PSCL0001	225	9	2,025

EA = each; OB = outer bag; PCH = pouch; CS = case

Notes

- a) The data shown are typical values and should not be used as product specifications.
b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.
d) All of Contec's packaging is compatible with hydrogen peroxide gassing applications.
e) These wipes are free of lint and loose fibers, and meet the definition of lint-free/low linting wipes according to the United States Pharmacopoeia Chapter 797 (USP-NF General Chapter <797> Pharmaceutical Compounding - Sterile Preparations) and the Institute of Environmental Sciences and Technology Recommended Practice IEST-RP-CC004.4

Test Methods:

- CTM = Contec Test Method
- IEST-RP-CC004.3 = Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of Environmental Sciences and Technology, Rolling Meadows IL.

PDSW043 | 052005

Copyright © 2022 Contec, Inc. All rights reserved.

Contact us:
info@contecinc.com
contecinc.com

Contec, Inc.
Spartanburg, SC 29303 USA
tel: +1-864-503-8333

