

Contec® Sterile IPA

Sterile isopropanol in purified water

Contec® Sterile isopropyl alcohol is a ready-to-use IPA solution pre-filtered to 0.2 micron and gamma irradiated. Containing USP Grade IPA and USP grade purified water and is designed for use in pharmaceutical, medical device, and other life science industries. Each bottle is double-bagged to ensure sterility and printed with lot number and date of expiration for easy record keeping. Validated sterile to a 10⁻⁶ SAL per ANSI/AAMI/ISO 11137 guidelines.

Note: Sterility validation is performed after trigger sprayer is installed.

A convenient trigger sprayer with settings for stream or mist spraying is included with 16 and 32 oz. bottles. The stream setting is ideal for targeting specific hard to reach areas, while the mist setting is used to cover larger areas. Each leakproof bottle is induction sealed, capped and shipped with an uninstalled trigger sprayer to eliminate the occurrence of leakage. Instructions for sterile installation of the trigger sprayer are available.



Product Application

- Wipe down of isolators, cabinets, laminar flow hoods and surfaces
- Use on gloved hands in sterile suites
- Wipe down products and equipment for pass through to controlled environments
- Ideally suited for use with sterile cleanroom wipe

Features	Benefits
Contains USP grade isopropyl alcohol	• Evaporates and leaves extremely low residues • Eliminates mixing, filtration, processing and QC expense of in-house production of sterile isopropyl alcohol water blends • Compliant with USP <797> and USP <800>
Filtered through 0.2 micron in-line filter	• Eliminates the need for in-house filtration
Gamma irradiated to a sterility assurance level (SAL) of 10 ⁻⁶	• Prefiltered, packaged and sterilized for easy introduction into critical environments • Certificate of processing and sterility with each case
Plastic bottle	• Recyclable, environmentally friendly, no costly disposal

Part No.	Description		Size	Packaging
SB167030IR	Contec Sterile 70% IPA Sterile 70% Isopropyl alcohol with trigger sprayer		16 oz. (0.47 L)	12/case
SB327030IR	Contec Sterile 70% IPA Sterile 70% Isopropyl alcohol with trigger sprayer		32 oz. (0.95 L)	12/case
SB1287030IR	Contec Sterile 70% IPA Sterile 70% Isopropyl alcohol		1 gallon (3.8 L)	4/case
SB128100IR	Contec Sterile IPA Sterile 99% Isopropyl alcohol		1 gallon (3.8 L)	4/case

Product Information		
Appearance	Clear, colorless liquid	
Odor	Characteristic alcohol	
Water Quality	Meets the requirements of USP grade deionized water	
% IPA Contents		
70%	68.0% - 72.0%	
100%	99.0% - 99.99%	
Specific Gravity (20C)		
70%	0.872% - 0.883%	
100%	0.785% - 0.787%	
Packaging materials	Outer bags (OB1, OB2), low density polyethylene (LDPE)  Case (CS), corrugated fiberboard (PAP) 	
Environment	ISO 3-8 Grade A/B	
Shelf life		
Unopened	5 years (SB167030IR, SB327030IR, SB1287030IR, SB128100IR)	
Opened	6 months (all)	

Recycle Symbols

PET	
HDPE	
LDPE	
PP	
PAP	

VOC Content	VOC (LBS/CS)	VOC (LBS/PCH)
SB167030IR	0.72	0.06
SB327030IR	1.51	0.12
SB1287030IR	1.99	0.49
SB128100IR	25.99	6.49

Contec® Disinfectant and Solution Rotation

DAILY	Use everyday on all surfaces.	Disinfectants CyQuanol™ Sterile CyQuanol™ PREempt® Plus PREempt® Concentrate	
ROUTINELY	Use based on room environment monitoring.	Sporicide PeridoxRTU® Sterile PeridoxRTU Peridox® Concentrate	
PERIODICALLY	Use to remove built up residue or when needed.	Detergent NeutraKlean Sterile NeutraKlean	
RINSE			
	Use after cleaning and disinfecting.	Rinse IPA HydroKlean WFI Water	